



Biplanar® - Designed, Developed and Manufactured in Sweden



Biplanar[®] - Designed, Developed and Manufactured in Sweden



Swemacimaging



The Biplanar G-Arm



The Method



Sales

Session Expectations

- Good understanding of the Biplanar Product
- Good understanding of the Biplanar Method
- Good understanding of the Main Selling Points
- Get sales going



Swemacimaging

Our History...

SAAB Aerospace

- 50 years ago, Swemac together with a SAAB Aerospace incubator developed the first G-arm, the Multiplane™
- It was conceived by two surgeons in Linköping, Sweden with the help of SAABs renowned engineers



...Our Present

Scanflex Healthcare Group

- Part of the Scanflex Healthcare Group since 2003
- Continued development and growth for the future



Scanflex | *Medical*
Scanflex | *Sirad*



 **Swemacimaging™**



Scanflex | *Dematech*

Product Evolution




1958-Multiplane



2021-Biplanar® 600s

Markets

- 900+ Units installed
- Key Markets, Nordics, China & Japan
- Focus on expansion



Organizational structure


Our organization



International Distribution




Support given by Swemacimaging




The Biplanar

G-Arm

The Biplanar Stand




The Biplanar Control Unit



The Biplanar Control Monitor




<table><tr><td>Patient</td><td>Patient-ID</td></tr><tr><td>Test Christa</td><td>37</td></tr><tr><td>Michaelis Heidrun</td><td>208</td></tr><tr><td>Orthmann Heinz</td><td>16</td></tr><tr><td>Allmann Gustav</td><td>147</td></tr></table>		Patient	Patient-ID	Test Christa	37	Michaelis Heidrun	208	Orthmann Heinz	16	Allmann Gustav	147	Hospital/Department																																
Patient	Patient-ID																																											
Test Christa	37																																											
Michaelis Heidrun	208																																											
Orthmann Heinz	16																																											
Allmann Gustav	147																																											
		Surgeon																																										
		Patient Name																																										
		Patient ID																																										
		Date of Birth																																										
		Accession Number																																										
		Image Description																																										
		System Messages																																										
<table><tr><td>kVp</td><td>mA</td><td>Pat.dose mGycm²</td><td>X-ray time</td><td>X-ray view</td><td>Noise suppression</td></tr><tr><td></td><td></td><td></td><td></td><td>☒ Blue=AP</td><td>☐ 1</td></tr><tr><td></td><td></td><td></td><td></td><td>☐ Yellow=lateral</td><td>☐ 2</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 4</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 6</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 8</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 16</td></tr></table>		kVp	mA	Pat.dose mGycm²	X-ray time	X-ray view	Noise suppression					☒ Blue=AP	☐ 1					☐ Yellow=lateral	☐ 2						☐ 4						☐ 6						☐ 8						☐ 16	Dicom
kVp	mA	Pat.dose mGycm²	X-ray time	X-ray view	Noise suppression																																							
				☒ Blue=AP	☐ 1																																							
				☐ Yellow=lateral	☐ 2																																							
					☐ 4																																							
					☐ 6																																							
					☐ 8																																							
					☐ 16																																							
☐ Contrast enhancement		☐ Edge enhancement		☐ Pos/Neg																																								
New Patient (F1)	Image Save (F2)	Image search (F3)	Image Delete (F4)	Pos/Neg (F5)	MediaPlayer (Ctrl+F6)	- Noise suppr. (F7)	+ Noise suppr. (F8)	User Man (Ctrl+F9)																																				

		Hospital/Department																																										
		Surgeon																																										
		Patient Name																																										
		Michaelis Heidrun																																										
		Patient ID																																										
		208																																										
		Date of Birth																																										
		12-08-1969																																										
		Accession Number																																										
		20																																										
		Image Description																																										
		System Messages																																										
<table><tr><td>kVp</td><td>mA</td><td>Pat.dose mGycm²</td><td>X-ray time</td><td>X-ray view</td><td>Noise suppression</td></tr><tr><td>85</td><td>2.3</td><td>00006</td><td>2.05</td><td>☒ Blue=AP</td><td>☐ 1</td></tr><tr><td></td><td></td><td></td><td></td><td>☐ Yellow=lateral</td><td>☐ 2</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 4</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 6</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 8</td></tr><tr><td></td><td></td><td></td><td></td><td></td><td>☐ 16</td></tr></table>		kVp	mA	Pat.dose mGycm²	X-ray time	X-ray view	Noise suppression	85	2.3	00006	2.05	☒ Blue=AP	☐ 1					☐ Yellow=lateral	☐ 2						☐ 4						☐ 6						☐ 8						☐ 16	Dicom
kVp	mA	Pat.dose mGycm²	X-ray time	X-ray view	Noise suppression																																							
85	2.3	00006	2.05	☒ Blue=AP	☐ 1																																							
				☐ Yellow=lateral	☐ 2																																							
					☐ 4																																							
					☐ 6																																							
					☐ 8																																							
					☐ 16																																							
☐ Contrast enhancement		☐ Edge enhancement		☐ Pos/Neg																																								
New Patient (F1)	Image Save (F2)	Image search (F3)	Image Delete (F4)	Pos/Neg (F5)	MediaPlayer (Ctrl+F6)	- Noise suppr. (F7)	+ Noise suppr. (F8)	User Man (Ctrl+F9)																																				

The Biplanar Control Panel





AP image settings




Lateral image settings

Common settings


Remote control and Foot switch



Options – What are they



MPPS - Modality Performed Procedure Step



RIS – Radiology Information System

PACS – Picture Archive and Communication System

The Method

Application & Method

Applications

- Hip surgery
- Femur surgery
- Tibia surgery
- Spine surgery
 - Kyphoplasty / vertebroplasty
 - Artificial disc replacement
 - Insertion of pedicle screws
 - Pain management
- Extremities surgery (ankle and foot)



The Biplanar Method Advantages

- Rotation free
- Images always in a 90 degree angle
- Stationary during surgery
- Patient is rotated if needed
- Minimize unnecessary fluoroscopy
- Speed & Accuracy
- Free up time for support staff




Positioning





Movements of the Biplanar™

- In the horizontal plane
 - pushing it around- light and runs on wheels
 - Each side has a pivot brake → rotation around a central point
- In the vertical plane → rotation (15 degrees in each direction) and elevation (25 cm)

Table

- All fluoroscopy tables can be used with the Biplanar
 - Certain tables like Jackson block the Biplanar from rolling under it
- If necessary, the table can also be moved





Examples

Intratrochanteric fracture –
Twin hook

Biplanar 500e



Examples

Kyphoplasty 2018

Biplanar 500e



Sales


Sales Focus

We are selling a tool that supports a method

- ❖ Hospital KPI's and targeted procedures
 - Dose, time, quality, recovery...
 - Orthopedic / Trauma / Pain management departments
- ❖ University hospitals
 - Highest budget available
 - Highest throughput of residents, spread the Biplanar® method
- ❖ Spine Center: large budget, great focus

- ❖ Find a Key Opinion Leader: test the Biplanar®, spread the method, write articles, make a video...
- ❖ Get a reference installation: bring doctors from other hospitals, participate in/attend surgeries

Pricing is secondary to the value



Sales arguments

Advantages

- ❖ Rotation-free AP and lateral imaging
- ❖ Simultaneous AP and lateral imaging
- ❖ Shorter operating time, safer surgery
- ❖ Ease of use
- ❖ Mobile



Benefits

- ❖ High degree of precision during surgery
- ❖ Less radiation
- ❖ Less risk of compromising the sterility
- ❖ Less anaesthesia
- ❖ Less bleeding
- ❖ Safer for the patient
- ❖ Better working conditions for the OR-staff



Return on Investment (ROI) Calculation

Assumptions

1 - At least 20 images are taken during a surgery (low figure)

2 - Time saved with Biplanar® is approx. 1 min / image (experience and feedback from surgeons)

3 - Procedure cost 20-100 € / min (varies from place to place)

4 - Five (5) procedures per week (many hospitals do 3 per day)

5 - 50 working weeks / year

Time & Quality = Savings


Saving of **20 min** for a normal procedure but up to 25% saving has been reported

Saving of **400-2 000 € /procedure**

Saving of **2 000–10 000 € / week**

Saving of **100 000–500 000 € / year**

We generally expect that a hospital that shifts to the Biplanar method will have a payback within **1-2 years**



Competition...

Competition cannot be viewed as any imaging device

Non-mobile imaging device for diagnostics...

Mobile imaging device, O-arm...

Non-mobile biplane imaging device for diagnostics...

Mobile imaging device for general surgery...

Mobile imaging device 3D C-arm...



Quiz

Quiz

- What are the main applications for the Biplanar?
- Name three advantages the Biplanar.
- What is the estimated yearly saving for a hospital implementing the Biplanar compared to a C-arm?



Biplanar[®] - Designed, Developed and Manufactured in Sweden